



AESONE CONNECT

DASHBOARD

- ✓ Accessible from your laptop, desktop & tablet devices
- ✓ Monitor & control your chamber in real-time
- ✓ Customizable features


SECURITY

- ✓ Control user access
- ✓ 3 authentication methods
- ✓ 6 permission levels



PROFILE

- ✓ Manage up to 10 profiles on your chamber
- ✓ Import and export profiles in JSON format



DATA LOGGING

- ✓ Monitor live & historical data
- ✓ Store up to 8 years of historical sensor data (dependent on how many sensors and what rate they are being logged)
- ✓ Pause & resume live plotting
- ✓ View, manage and export activity logs in JSON format
- Export sensor data in csv/tsv format

FIND OUT MORE:

www.AssociatedEnvironmentalSystems.com
sales@associatedenvironmentalsystems.com | 978.772.0022