

---

# Job Parts List

**Job #** BHD-405

---

<b>Part</b>	<b>Description</b>	<b>Qty Per</b>	<b>Pick / Buy</b>
BOX-MEDIUM	62" X 48" X 55"H	1	P
CAP01839-7	.042 Cap Tubing, 7'	1	P
CAP01842-6	.055 Cap Tubing X 6'	1	P
CAP03791-9	.049 Cap Tubing X 9'	1	P
CAP03999	Strainer	1	P
CAS02140	Leveling Glide Legs	0.5	P
CIB03158	Circuit Breaker 25A 125/250VAC	2	P
CIB03247	Circuit Breaker 5 Amp, 250 VAC	2	P
CIB03248	Circuit Breaker 15 Amp, 250VAC	4	P
CIB03250	Circuit Breaker 8 Amp, 250 VAC	2	P
CMP02530	Compressor CS14K6E-PFV-256	1	P
CMP02679	Compressor Mounting Kit	1	P
CMP02910	Capacitor Start	1	P
CMP02911	Capacitor Run	1	P
CMP02912	Relay Kit	1	P
CMP03451	Retaining Rings - Small	1	P
CMP03452	Retaining Rings - Large	1	P
CTL01502	Sensor, 100 OHM Platinum	1	P
CTL01828	Limit Switch, Automatic	1	P
CTL01829	Limit Switch, Manual	1	P
CTL02160-1	Type J Thermocouple, 72"	1	P

<b>Part</b>	<b>Description</b>	<b>Qty Per</b>	<b>Pick / Buy</b>
CTL02434	Sensor, Humidity 4-20mA	1	P
CTL02519	Controller - Two Channel	1	P
CTL03403	High Temperature Failsafe	1	P
DRY01739	Filter Drier	1	P
ELD03286	Quencharc Resistor, Paktron	2	P
ELD04224	Ground Bars - 10	1	P
ELE01673	Female Connector, 9 Pin	1	P
ELF02958	Strain Relief for 3/4" Cord	1	P
ELS01209	Switch 2P1T	6	P
ELS01228	Clear Cap	2	P
ELS01230	Amber Cap	4	P
ELW02009	Line Cord, 30A 220V	1	P
FBW01014	Fan Blade 6" Diameter	1	P
FBW02662	Fan Blade Condenser, 12"	1	P
FCL01913	Evaporator 10X19.25 1HP	1	P
FCL02663	Condenser 12 X 12.25	1	P
FSE01342	Heat Fuse, WUE	2	P
FSE02604	Fuse 1A Glass	3	P
FSE03186	Glass Fuse .75A	2	P
GSK04066	Grey Silicone Sponge Gasket	8	P
HDW01110	Handle, Door	1	P
HDW01118	4 X 4 Riveted Pin Hinge	2	P
HDW01129	Door Latch - Large	1	P
HDW01141	Door Strike - Large	1	P

<b>Part</b>	<b>Description</b>	<b>Qty Per</b>	<b>Pick / Buy</b>
HTR01245	Ceramic Heater - 1000W 240V	2	P
HTR01247	Heater, Humidity 600W 240V	1	P
HUM01540	Float Valve Assembly	1	P
HUM04137	Steam Generator, SS, BHD-8/5	1	P
LGT01231	Light Fixture	1	P
LGT03018	Light Bulb, 60W 130V Tubular	1	P
MTR02661	Motor Condenser	1	P
MTR03768	Circulation Motor Bracket 1/2"	1	P
MTR04248	Motor, Circ. 6-5/8" Shaft	1	P
PCB03226	Circuit Board, 220V	1	P
PRB01582	Port Sleeve, 1/2"	1	P
PRB01583	Port Sleeve, 7/8"	1	P
PRB01620	1/2" Port Collar	1	P
PRB01621	7/8" Port Collar	1	P
PRB03991	Silicone 3" Port	1	P
PRB04062	3" Sponge Plug	1	P
PRS03910	Pressure Switch, NC, 400 PSI	1	P
PTF03426	Tee, Nylon 3/8X3/8	1	P
PTF03427	Ball Valve, 2 Way PVC 3/8"	1	P
PTF03428	Adapter, Nylon, Female 3/8X1/2	1	P
REL02832	Off Delay Timer	1	P
REL03256	Relay, Crydom	4	P
SLV03096	Shelf Rails - 4 foot lengths	2	P
TRS01447	Transformer, 500VA 220V	1	P

<b>Part</b>	<b>Description</b>	<b>Qty Per</b>	<b>Pick / Buy</b>
TRS03258	Transformer, 24VAC	1	P
TXV02233	Valve, Expansion 3X4 ODF 30"	1	P
VAL01719	Solenoid Valve	3	P
VAL01721	Solenoid Coil 120V	3	P
VAL03906	Access Valve, 1/4"	2	P
WIN01168	Window 9" X 9"	1	P
WRH03272	Wire Harness J3-4	1	P
WRH03275	Wire Harness J17-8	1	P
WRH03277	Wire Harness J2-4	1	P
WRH03281	Wire Harness J23-8	1	P
WRH03284	Wire Harness J25-4	1	P
WRH03292	Wire Harness J8-8	1	P
WRH03294	Wire Harness J20-8	1	P
WRH03295	Wire Harness J21-8	1	P
WRH03297	Wire Harness J22-8	1	P
WRH03382	Wire Harness J10-8	1	P