
Job Parts List

Job # HD-464

Part	Description	Qty Per	Pick / Buy
ACU01750	Suction Accumulator, 1-1/8"	1	B
CAP03999	Strainer	1	P
CAS02150	Brake Kit	2	B
CAS02557	Caster, 5" Swivel	2	B
CAS02558	Caster, 5" Rigid	2	B
CMP01663	Sentronic Oil Safety	1	B
CMP02469	Compressor	1	B
CTL01502	Sensor, 100 OHM Platinum	1	P
CTL01828	Limit Switch Manual	3	P
CTL01829	Limit Switch, Automatic	3	P
CTL02160	Thermocouple Type J, 72"	1	P
CTL02434	Sensor, Humidity 4-20mA	1	P
CTL02519	Controller - Two Channel	1	P
CTL03403	High Temperature Failsafe	1	P
DIS01804	Auxiliary Connector	1	B
DRY01737	Filter Drier, 3/8	1	B
ELB04438	8 X 8 electrical box	1	P
ELD01799	Distribution Block	1	P
ELD03286	Quencharc Resistor, Paktron	1	P
ELE01670	25 Pin Female RS232 Connector	1	P
ELS01209	Switch 2P1T	6	P

Part	Description	Qty Per	Pick / Buy
ELS01228	Clear Cap	2	P
ELS01230	Amber Cap	4	P
FBW01041	Blower Wheel, SS, 7-11/16"	2	P
FBW02444	Fan Blade 14" 3 Blade	2	P
FCL02471	Condenser, 3HP Air Cooled	1	P
FCL02628	Evaporator 5 HP, B30-5697	1	P
FSE01268	2P 30A Fuse Block	14	P
FSE01273	Fuse 1A OT 250V	1	P
FSE01279	Fuse, 20A OT 250V	6	P
FSE01280	Fuse, 25A OT 250V	2	P
FSE01284	Fuse 1A Form 101	2	P
FSE01310	Fuse, 4A TD 250V	2	P
FSE01312	Fuse 6-1/4A TD 250V	3	P
FSE01317	Fuse 30A TD 250V	3	P
FSE01342	Heat Fuse	6	P
FSE01351	Fuse, 3.5A TD 250V	2	P
FSE02604	Fuse 1A Glass	2	P
FSE02605	Fuse Holder - Inline	2	P
GSK04066	Grey Silicone Sponge Gasket	16	P
HDW01103	Door Handle, Large	1	P
HDW01104	Latch, Black	1	P
HDW01148	6 X 6 Blank Hinge	0	P
HDW01149	Latch, Adjustable Grip	12	P
HTR01245	Ceramic Heater - 1000W 240V	12	P

Part	Description	Qty Per	Pick / Buy
HTR01248	Heater, Humidity 1500W 240V	3	P
HTR01253	Heater, Crankcase	1	B
HUM01540	Reservoir Assembly	1	P
HUM04138	Steam Generator, SS, HD-64	1	P
LGT02155	Light Fixture, Vapor Proof	1	P
LGT04068	75W Light Bulb	1	P
MTR01002	Motor, 1/3HP, 1PH, 50/60HZ	2	P
MTR02443	Motor, Condenser Fan 208V 1.6A	2	B
PRB01582	Port Sleeve, 1/2"	1	P
PRB01583	Port Sleeve, 7/8"	2	P
PRB01620	1/2" Port Collar	1	P
PRB01621	7/8" Port Collar	2	P
PRB03991	Silicone 3" Port	1	P
PRB04062	3" Sponge Plug	1	P
PRS02332	Pressure Control, Dual	1	B
REC01971	Receiver	1	B
REL01397	Relay 2P2T	2	P
REL02382	Solid State Relay 1P1T	4	P
REL02832	Off Delay Timer	1	P
REL03944	Blisterpak for Solid State	1	P
SEP02122	Oil Separator	1	B
SGL01184	Sight Glass, 3/8"	1	B
TRS01447	Transformer, 500VA 220V	1	P
TRS02373	Transformer, 10 VA	1	P

Part	Description	Qty Per	Pick / Buy
TXV01788	Valve, Expansion	1	B
TXV02719	Expansion Valve	1	B
VAL01552	Shut Off Valve	1	B
VAL01566	Solenoid Valve	2	B
VAL01706	Pressure Relief Valve	1	B
VAL01719	Solenoid Valve	1	P
VAL01721	Solenoid Coil 120V	1	P
VAL01722	Coil, Solenoid 120V	2	B
VAL01819	Valve, Rotalock	1	B
VAL01958	Hot Gas By Pass Valve	1	B
VAL02813	Valve, Packless Shut Off	1	B
VAL03906	Access Valve, 1/4"	2	P
VIB02127	Vibration Eliminator 7/8"	1	B
VIB02130	Vibration Eliminator, 1-1/8"	1	B
WIN01170	Window 13" X 13"	1	P