

---

# Job Parts List

**Job #** LH-6

---

<b>Part</b>	<b>Description</b>	<b>Qty Per</b>	<b>Pick / Buy</b>
CAP01837-8	.031 Cap Tubing, 8'	1	P
CAP01839-12	.042 Cap Tubing, 12'	1	P
CAP01842-4	.055 Cap Tubing, 4'	2	P
CAP03999	Strainer	1	P
CAS02140	Leveling Glide Legs	0.5	P
CIB03158	Circuit Breaker 25A 125/250VAC	1	P
CIB03248	Circuit Breaker 15 Amp, 250VAC	2	P
CIB03250	Circuit Breaker 8 Amp, 250 VAC	2	P
CMP02281	ECL Condensing Unit	1	P
CTL01502	Sensor, 100 OHM Platinum	1	P
CTL01828	Limit Switch Manual	1	P
CTL01829	Limit Switch, Automatic	1	P
CTL02161	Failsafe Control	1	P
CTL03338	PM Controller	2	P
CTL03378	RH Sensor for LH-1.5	1	P
DRY01741	Filter Drier	1	P
ELC03377	Receptacle, Duplex 15A 125V	1	P
ELC03726	Power Supply 24VDC	1	P
ELD01363	Jumper	2	P
ELD03286	Quencharc Resistor, Paktron	2	P
ELF02868	Strain Relief for 14-3 Line	0	P

<b>Part</b>	<b>Description</b>	<b>Qty Per</b>	<b>Pick / Buy</b>
ELS01209	Switch 2P1T	2	P
ELS01228	Clear Cap	1	P
ELS01230	Amber Cap	1	P
ELW02008	Line Cord 15A 125V 6' 14/3 SO	0	P
FBW01014	Fan Blade 6" Diameter	1	P
FCL01912	Evaporator 1/3 HP, 8 X 11	1	P
FSE02604	Fuse 1A Glass	1	P
FSE03186	Glass Fuse .75A	2	P
GSK04066	Grey Silicone Sponge Gasket	8	P
HDW01106	Latch for ECL	1	P
HDW01115	Hinge, 3-1/2 X 3-1/2"	2	P
HTR01251	Heater, Humidity 600W 120V	1	P
HTR03936	Ceramic Heaters, 500W 120V	2	P
HUM01540	Reservoir Assembly	1	P
MTR01001	Motor 6-5/8 Shaft	1	P
PCB03227	Circuit Board, 110V	1	P
PRB01582	Port Sleeve, 1/2"	1	P
PRB01583	Port Sleeve, 7/8"	2	P
PRB01620	1/2" Port Collar	1	P
PRB01621	7/8" Port Collar	2	P
PRB03991	Silicone 3" Port	1	P
PRB04062	3" Sponge Plug	1	P
PRS03910	Pressure Control for LH Series	1	P
PTF03426	Tee, Nylon 3/8X3/8	1	P

<b>Part</b>	<b>Description</b>	<b>Qty Per</b>	<b>Pick / Buy</b>
PTF03427	Ball Valve, 2 Way PVC 3/8"	1	P
PTF03428	Adapter, Nylon, Female 3/8X1/2	1	P
REL03256	Relay, Crydom	2	P
SLV02899	Shelf Clips, Stainless Steel	4	P
SLV03096	Shelf Rails	2	P
SLV03994	Shelf, 17 X 19-1/2" for LH-6	1	P
TRS03793	24V Transformer Bracket	1	P
VAL01719	Solenoid Valve	2	P
VAL01721	Solenoid Coil 120V	2	P
VAL03906	Access Valve, 1/4"	2	P
WRH03272	Wire Harness J3-4	1	P
WRH03277	Wire Harness J2-4	1	P
WRH03278	Wire Harness J20-4	1	P
WRH03281	Wire Harness J23-8	1	P
WRH03292	Wire Harness J8-8	1	P
WRH03295	Wire Harness J21-8	1	P
WRH03297	Wire Harness J22-8	1	P