
Job Parts List

Job # MX-9208

Part	Description	Qty Per	Pick / Buy
BOX-MEDIUM	62" X 48" X 55"H	1	P
ELF02868	Strain Relief for 14-3 Line	1	P
ELS01522	Toggle Switch, 1P1T	3	P
ELS02422	Pilot Light - Circulation	1	P
ELS02925	Pilot Light - Chamber	1	P
ELS02926	Pilot Light - Tower	1	P
ELW02008	Line Cord 15A 125V 6' 14/3 SO	1	P
FBW01023	Blower, Salt Spray	1	P
FSE01268	Fuse Block, 2P 30A 250V	1	P
FSE01278	Fuse 15A OT 250V	1	P
HDW04107	1/4" x 1/4" Nylon Fitting	14	P
HDW04108	1/4" x 1/4" Nylon Elbow	5	P
HDW04109	1/4" X 1/4" Nylon Bulkhead	1	P
HDW04117	1/4" Nylon Plug	3	P
HDW04118	1/4" x 1/8" Nylon Elbow	1	P
PRB01609	Recess Bumper	4	P
PRB01628	2-1/8" Locking Ring	1	P
SSC02563	Heater, 250W Immersion	1	P
SSC02564	Thermometer, Salt Spray, Tower	1	P
SSC02565	Thermostat	2	P
SSC02566	2-1/4" Super Vacuflex Hose	6	P

Part	Description	Qty Per	Pick / Buy
SSC02567	Level Switch, Salt Spray	1	P
SSC02569	Check Valve 1/4"NPT 1 PSI	1	P
SSC02575	Gasket: Tower	2	P
SSC02576	Heater, 500W Tubular Air	1	P
SSC02577	Valve, Ball 1/4" FPT 2 Way PVC	4	P
SSC02578	Thermometer, Salt Spray, Cover	1	P
SSC02579	Valve, 3 Way Brass	1	P
SSC02580	Relief Valve, Salt Spray	1	P
SSC02581	Pressure Gauge 2", Salt Spray	1	P
SSC02582	Clamp: SS, 2.06-2.75"	6	P
SSC02873	Siphon Tube 1/2" X 18 ³ / ₄ "	1	P
SSC02883	Solenoid Valve	1	P
SSC02886	3/4" Slip PVC Elbow	1	P
SSC02887	3/4" Male Adapters	1	P
SSC02889	Reducing Coupling DMV	1	P
SSC02927	Pressure Switch, N/O-C	1	P
SSC02928	Tower Floor for MX-9204	1	P
SSC02930	Tower Top	1	P
SSC02931	PVC Pipe 3/4", Schedule 40	1	P
SSC02932	Tower Rod - Length 21"	4	P
SSC02934	Nut, Acorn 1/4-20 SS	4	P
SSC02935	Filter: Nylon Cloth	1	P
SSC02936	Grommet: Rubber	2	P
SSC02975	Salt Spray Lucite Assembly	1	P

Part	Description	Qty Per	Pick / Buy
SSC03835	Plastic Pipe Nipple	1	P
SSC03836	Plastic Pipe Elbow	1	P
SSC03838	Pressure Switch	1	P
SSC03851	1/8 X Close SS 304	1	P
SSC03852	1/8" x 1/8" SS 304 Coupling	1	P
SSC03902-1	Salt Spray Nozzle Body	1	P
SSC03902-2	Salt Spray Nozzle, Large 8+	1	P
SSC03902-3	Salt Spray Nozzle, Ret. Ring	1	P
SSC03902-4	Salt Spray Nozzle, Gasket	1	P
SSC03921	Hanging Lucite Rod	6	P
SSC03924	Lucite Tank Plate	1	P