
Job Parts List

Job # BD-100

Part	Description	Qty Per	Pick / Buy
CTL01502	Sensor, 100 OHM Platinum	1	P
CTL02160	Thermocouple Type J, 48"	1	P
CTL03338	PM Controller	1	P
CTL03403	High Temperature Failsafe	1	P
ELD01362	Quencharc Resistor, Paktron	1	P
ELD03179	Miniblock, Terminal Board	4	P
ELD03180	Miniblock Terminal Block	2	P
ELD03181	Miniblock End Section,	1	P
ELD03182	Miniblock End Cap Section	1	P
ELE01673	Female Connector, 9 Pin	1	P
ELF02868	Strain Relief for 14-3 Line	1	P
ELS01209	Switch 2P1T	3	P
ELS01228	Clear Cap	1	P
ELS01230	Amber Cap	2	P
ELW02008	Line Cord 15A 125V 6' 14/3 SO	1	P
FBW01013	4" Fan Blade	1	P
FIL02545	CO-2 Scrubber Housing	1	P
FIL02546	CO-2 Scrubber Cartridge -6 pac	0.17	P
FSE01268	2P 30A Fuse Block	2	P
FSE01278	Fuse 15A OT 250V	1	P
FSE01284	Fuse 1A Form 101	1	P

Part	Description	Qty Per	Pick / Buy
FSE01311	Fuse, 5A, TD, 250V	1	P
FSE01314	15A TD 250V Fuses	1	P
FSE02604	Fuse 1A Glass	2	P
FSE02605	Fuse Holder - Inline	2	P
FSE03177	Heat Fuse, High Temp.	1	P
GSK01157	3/8" X 1" Grey Silicone Sponge	5	P
HDW01110	Handle, Door	1	P
HDW01115	Hinge, 3-1/2 X 3-1/2"	2	P
HDW01128	Latch, Door Small	1	P
HTR01242	660W 120V Ceramic Heater	1	P
MTR01001	Motor 6-5/8 Shaft	1	P
PRB01582	Port Sleeve, 1/2"	3	P
PRB01583	Port Sleeve, 7/8"	2	P
PRB01598	Port Sleeve 2" OD	1	P
PRB01599	2" Port Plug	1	P
PRB01609	Recess Bumper	4	P
PRB01620	1/2" Port Collar	3	P
PRB01621	7/8" Port Collar	2	P
PRB01627	Locking Ring 2"	1	P
REL01408	Contactora 2P 30A 120V Coil	1	P
REL02382	Solid State Relay 1P1T	1	P
REL03943	MOV For Solid State Relay	1	P
REL03944	Blisterpak for Solid State	1	P
VAL02028	Valve, CO-2 1000 PSI 120V	1	B

Part	Description	Qty Per	Pick / Buy
VAL02040	Orifice, CO-2 1000 PSI	1	P
VAL02052	Pipe Nipple, SS 1/8" X 5"L	1	P