
Job Parts List

Job # BD-108

Part	Description	Qty Per	Pick / Buy
CTL01502	Sensor, 100 OHM Platinum	1	P
CTL02160	Thermocouple Type J, 48"	1	P
CTL03338	PM Controller	1	P
CTL03403	High Temperature Failsafe	1	P
ELF02958	Strain Relief for 3/4" Cord	1	P
ELS01209	Switch 2P1T	3	P
ELS01228	Clear Cap	1	P
ELS01230	Amber Cap	2	P
ELW02011	Line Cord, 10/3 SJT L5-30P	1	P
FBW01014	Fan Blade 6" Diameter	1	P
FIL02545	CO-2 Scrubber Housing	1	P
FIL02546	CO-2 Scrubber Cartridge -6 pac	0.127	P
FSE01268	2P 30A Fuse Block	2	P
FSE01279	Fuse, 20A OT 250V	1	P
FSE01284	Fuse 1A Form 101	1	P
FSE01311	Fuse, 5A, TD, 250V	1	P
FSE01316	Fuse, 25A 250V TD	1	P
FSE02604	Fuse 1A Glass	2	P
FSE02605	Fuse Holder - Inline	2	P
FSE03177	Heat Fuse	1	B
GSK04066	Grey Silicone Sponge Gasket	10	P

Part	Description	Qty Per	Pick / Buy
HDW01110	Handle, Door	1	P
HDW01118	4 X 4 Riveted Pin Hinge	2	P
HDW04181	LRG Adj Draw Pull Latch, SS	1	P
HDW04182	Keeper/Strike for 2173, 304SS	1	P
HTR01242	660W 120V Ceramic Heater	3	P
MTR01001	Motor 6-5/8 Shaft	1	P
PRB01582	Port Sleeve, 1/2"	3	P
PRB01583	Port Sleeve, 7/8"	2	P
PRB01620	1/2" Port Collar	3	P
PRB01621	7/8" Port Collar	2	P
PRB03991	Silicone 3" Port	1	P
PRB04062	3" Sponge Plug	1	P
REL01408	Contactactor 2P 30A 120V Coil	1	P
REL02382	Solid State Relay 1P1T	1	P
REL03943	MOV For Solid State Relay	1	P
REL03944	Blisterpak for Solid State	1	P
VAL02028	Valve, CO-2 1000 PSI 120V	1	B
VAL02040	Orifice, CO-2 1000 PSI	1	P
VAL02996	Pipe Nipple, SS 1/8" X 2" L	1	P