
Job Parts List

Job # BD-127

| Part | Description | Qty Per | Pick / Buy |
|-------------|--------------------------------|----------------|-------------------|
| CTL01502 | Sensor, 100 OHM Platinum | 1 | P |
| CTL03110 | Type T Thermocouple | 1 | P |
| CTL03338 | PM Controller | 1 | P |
| CTL03403 | High Temperature Failsafe | 1 | P |
| ELD01800 | Distribution Block, 2P | 1 | P |
| ELD03179 | Miniblock, Terminal Board | 6 | P |
| ELD03182 | Miniblock End Cap Section | 2 | P |
| ELE01673 | Female Connector, 9 Pin | 1 | P |
| ELS01209 | Switch 2P1T | 3 | P |
| ELS01228 | Clear Cap | 1 | P |
| ELS01230 | Amber Cap | 3 | P |
| FBW01041 | Blower Wheel, SS, 7-11/16" | 1 | P |
| FIL02545 | CO-2 Scrubber Housing | 1 | P |
| FIL02546 | CO-2 Scrubber Cartridge -6 pac | 0.17 | P |
| FSE01268 | 2P 30A Fuse Block | 5 | P |
| FSE01281 | Fuse, 30A OT 250V | 4 | P |
| FSE01284 | Fuse 1A Form 101 | 1 | P |
| FSE01312 | Fuse 6-1/4A TD 250V | 1 | P |
| FSE01351 | Fuse, 3.5A TD 250V | 4 | P |
| FSE02604 | Fuse 1A Glass | 2 | P |
| FSE02605 | Fuse Holder - Inline | 2 | P |

| Part | Description | Qty Per | Pick / Buy |
|-------------|-------------------------------|----------------|-------------------|
| FSE03177 | Heat Fuse, High Temp. | 4 | B |
| GSK04066 | Grey Silicone Sponge Gasket | 25 | P |
| HDW01110 | Handle, Door | 1 | P |
| HDW01122 | Hinge, Blank 5 X 5 | 2 | P |
| HDW04181 | LRG Adj Draw Pull Latch, SS | 1 | P |
| HDW04182 | Keeper/Strike for 2173, 304SS | 1 | P |
| HTR01245 | Ceramic Heater - 1000W 240V | 8 | P |
| MTR01002 | Motor, 1/3HP, 1PH, 50/60HZ | 1 | P |
| PRB01582 | Port Sleeve, 1/2" | 2 | P |
| PRB01583 | Port Sleeve, 7/8" | 1 | P |
| PRB01620 | 1/2" Port Collar | 2 | P |
| PRB01621 | 7/8" Port Collar | 1 | P |
| PRB03991 | Silcone 3" Port | 1 | P |
| PRB04062 | 3" Sponge Plug | 1 | P |
| REL01400 | Relay SSR, DC, 1P1T | 1 | P |
| REL01408 | Contactactor 2P 30A 120V Coil | 3 | P |
| REL02382 | Solid State Relay 1P1T | 4 | P |
| REL03943 | MOV For Solid State Relay | 4 | P |
| REL03944 | Blisterpak for Solid State | 1 | P |
| TRS01447 | Transformer, 500VA 220V | 1 | P |
| VAL02028 | Valve, CO-2 1000 PSI 120V | 1 | B |
| VAL02040 | Orifice, CO-2 1000 PSI | 1 | P |
| VAL02052 | Pipe Nipple, SS 1/8" X 5"L | 1 | P |