

---

# Job Parts List

**Job #** BHD-502

---

<b>Part</b>	<b>Description</b>	<b>Qty Per</b>	<b>Pick / Buy</b>
CAP01837-9	.031 Capillary Tubing X 9'	1	P
CAP01839-10	.042 Capillary Tubing X 10'	1	P
CAP01839-6	.042 Capillary Tubing X 6'	1	P
CAP03999	Strainer	2	P
CIB03247	Circuit Breaker 5 Amp, 250 VAC	2	P
CIB03248	Circuit Breaker 15 Amp, 250VAC	6	P
CIB03250	Circuit Breaker 8 Amp, 250 VAC	1	P
CMP02670	Compressor	1	P
CTL01502	Sensor, 100 OHM Platinum	2	P
CTL01828	Limit Switch, Automatic	1	P
CTL01829	Limit Switch, Manual	1	P
CTL02160-1	Type J Thermocouple, 72"	1	P
CTL02434	Sensor, Humidity 4-20mA	1	P
CTL02519	Controller - Two Channel	1	P
CTL03403	High Temperature Failsafe	1	P
DRY01739	Filter Drier	1	P
ELD03286	Quencharc Resistor, Paktron	2	P
ELF02958	Strain Relief for 3/4" Cord	1	P
ELS01209	Switch 2P1T	6	P
ELS01228	Clear Cap	2	P
ELS01230	Amber Cap	4	P

<b>Part</b>	<b>Description</b>	<b>Qty Per</b>	<b>Pick / Buy</b>
ELW02009	Line Cord, 30A 220V	1	P
FBW01014	Fan Blade 6" Diameter	1	P
FBW02867	Condenser Fan Blade 10"	1	P
FCL01912	Evaporator 1/3 HP, 8 X 11	1	P
FCL01926	Condenser 12 X 10.38 1 HP	1	P
FSE01342	Heat Fuse	2	P
FSE02604	Fuse 1A Glass	3	P
FSE03186	Glass Fuse .75A	2	P
GSK01157	3/8" X 1" Grey Silicone Sponge	7	P
HDW01110	Handle, Door	1	P
HDW01118	4 X 4 Riveted Pin Hinge	2	P
HDW01129	Door Latch - Large	1	P
HDW01141	Door Strike - Large	1	P
HTR01242	660W 120V Ceramic Heater	2	P
HTR01251	Heater, Humidity 600W 120V	1	P
HUM01540	Reservoir Assembly	1	P
LGT01231	Light Fixture	1	P
LGT03018	Light Bulb, 60W 130V Tubular	1	P
MTR01001	Motor 6-5/8 Shaft	1	P
MTR02661	Motor Condenser	1	P
MTR03327	Condenser Motor Bracket	1	P
PCB03227	Circuit Board, 110V	1	P
PRB01582	Port Sleeve, 1/2"	1	P
PRB01583	Port Sleeve, 7/8"	2	P

<b>Part</b>	<b>Description</b>	<b>Qty Per</b>	<b>Pick / Buy</b>
PRB01598	Port Sleeve 2" OD	1	P
PRB01599	2" Port Plug	1	P
PRB01620	1/2" Port Collar	1	P
PRB01621	7/8" Port Collar	2	P
PRB01627	Locking Ring 2"	1	P
PRS03910	Pressure Control, NC, 400 PSI	1	P
PTF03426	Tee, Nylon 3/8X3/8	1	P
PTF03427	Ball Valve, 2 Way PVC 3/8"	1	P
PTF03428	Adapter, Nylon, Female 3/8X1/2	1	P
REL02832	Off Delay Timer	1	P
REL03256	Relay, Crydom	2	P
SLV02899	Shelf Clips, Stainless Steel	4	P
SLV03096	Shelf Rails	4	P
TRS03258	Transformer, 24VAC	1	P
TXV02233	Valve, Expansion	1	P
VAL01719	Solenoid Valve	3	P
VAL01721	Solenoid Coil 120V	3	P
VAL03906	Access Valve, 1/4"	2	P
WIN01168	Window 9" X 9"	1	P
WRH03275	Wire Harness J17-8	1	P
WRH03281	Wire Harness J23-8	1	P
WRH03284	Wire Harness J25-4	1	P
WRH03289	Wire Harness J2-8	1	P
WRH03290	Wire Harness J3-8	1	P

<b>Part</b>	<b>Description</b>	<b>Qty Per</b>	<b>Pick / Buy</b>
WRH03292	Wire Harness J8-8	1	P
WRH03294	Wire Harness J20-8	1	P
WRH03295	Wire Harness J21-8	1	P
WRH03297	Wire Harness J22-8	1	P