
Job Parts List

Job # HD-405

Part	Description	Qty Per	Pick / Buy
CAP01854	Strainer/Screen	3	P
CAS02786	Caster Wheel, 3-1/2" Swivel	2	P
CAS02787	Caster Wheel, 3-1/2" Rigid	2	P
CAS02788	Caster Wheel Brake Kit	2	P
CMP02530	Compressor	1	P
CTL01502	Sensor, 100 OHM Platinum	1	P
CTL01828	Limit Switch, Automatic	1	P
CTL01829	Limit Switch, Manual	1	P
CTL02160	Thermocouple Type J, 48"	1	P
CTL02434	Sensor, Humidity 4-20mA	1	P
CTL02519	Controller - Two Channel	1	P
CTL03403	High Temperature Failsafe	1	P
DRY01739	Filter Drier	1	P
ELE01670	25 Pin Female RS232 Connector	1	P
ELF02958	Strain Relief for 3/4" Cord	1	P
ELS01209	Switch 2P1T	6	P
ELS01228	Clear Cap	2	P
ELS01230	Amber Cap	4	P
ELW02009	Line Cord, 30A 220V	1	P
FBW01014	Fan Blade 6" Diameter	1	P
FBW02662	Fan Blade Condenser, 12"	1	P

Part	Description	Qty Per	Pick / Buy
FCL01913	Evaporator 10X19.25 1HP	1	P
FCL02663	Condenser 12 X 12.25	1	P
FSE01268	2P 30A Fuse Block	7	P
FSE01273	Fuse 1A OT 250V	1	P
FSE01276	Fuse, 5A OT 250V	2	P
FSE01278	Fuse 15A OT 250V	2	P
FSE01284	Fuse 1A Form 101	2	P
FSE01312	Fuse 6-1/4A TD 250V	1	P
FSE01314	15A TD 250V Fuses	2	P
FSE01317	Fuse 30A TD 250V	2	P
FSE01342	Heat Fuse	2	P
FSE01351	Fuse, 3.5A TD 250V	2	P
FSE02604	Fuse 1A Glass	2	P
FSE02605	Fuse Holder - Inline	2	P
GSK01157	3/8" X 1" Grey Silicone Sponge	10	P
HDW01110	Handle, Door	1	P
HDW01118	4 X 4 Riveted Pin Hinge	2	P
HDW01129	Door Latch - Large	1	P
HDW01141	Door Strike - Large	1	P
HTR01245	Ceramic Heater - 1000W 240V	2	P
HTR01247	Heater, Humidity 600W 240V	1	P
HUM01540	Float Valve Assembly	1	P
LGT01231	Light Fixture	1	P
MTR01001	Motor 6-5/8 Shaft	1	P

Part	Description	Qty Per	Pick / Buy
MTR02661	Motor Condenser	1	P
PRB01582	Port Sleeve, 1/2"	1	P
PRB01583	Port Sleeve, 7/8"	2	P
PRB01598	Port Sleeve 2" OD	1	P
PRB01599	2" Port Plug	1	P
PRB01620	1/2" Port Collar	1	P
PRB01621	7/8" Port Collar	2	P
PRB01627	Locking Ring 2"	1	P
PRS01822	Pressure Switch	1	P
REL01397	Relay 2P2T	3	P
REL01408	Contactora 2P 30A 120V Coil	2	P
REL02382	Solid State Relay 1P1T	4	P
REL02832	Off Delay Timer	1	P
TRS01447	Transformer, 500VA 220V	1	P
TRS02373	Transformer, 24 VA	1	P
TXV02233	Valve, Expansion 3X4 ODF 30"	1	P
VAL01704	Access Valves	2	P
VAL01719	Solenoid Valve	3	P
VAL01721	Solenoid Coil 120V	3	P
WIN01168	Window 9" X 9"	1	P