
Job Parts List

Job # SD-604

Part	Description	Qty Per	Pick / Buy
CTL03235	High Temperature Sensor	1	P
CTL03338	PM Controller	1	P
CTL03340	PM High/Low Limit Control	1	B
ELD03286	Quencharc Resistor, Paktron	1	P
ELF02958	Strain Relief for 3/4" Cord	1	P
ELS01209	Switch 2P1T	3	P
ELS01228	Clear Cap	1	P
ELS01230	Amber Cap	2	P
ELW02009	Line Cord, 30A 220V	1	P
FBW01020	Axial Fan	1	P
FBW01022	Fan Blade, 6-1/2" SS	1	P
FBW03135	Fan Guard	1	P
FBW03203	Axial Fan Cordset	1	B
FIL02545	CO-2 Scrubber Housing	1	P
FIL02546	CO-2 Scrubber Cartridge -6 pac	0.17	P
FSE01268	2P 30A Fuse Block	4	P
FSE01279	Fuse, 20A OT 250V	2	P
FSE01284	Fuse 1A Form 101	1	P
FSE01312	Fuse 6-1/4A TD 250V	1	P
FSE01316	Fuse, 25A 250V TD	2	P
FSE01351	Fuse, 3.5A TD 250V	2	P

Part	Description	Qty Per	Pick / Buy
GSK01162	High Temperature Gasket	5	P
HDW01110	Handle, Door	1	P
HDW01123	Hinge, Riveted Pin 3" X 3"	2	P
HDW01129	Door Latch - Large	1	P
HDW01141	Door Strike - Large	1	P
HTR01245	Ceramic Heater - 1000W 240V	3	P
MTR01001	Motor 6-5/8 Shaft	1	P
PRB01582	Port Sleeve, 1/2"	2	P
PRB01583	Port Sleeve, 7/8"	2	P
PRB01598	Port Sleeve 2" OD	1	P
PRB01620	1/2" Port Collar	2	P
PRB01621	7/8" Port Collar	2	P
PRB01627	Locking Ring 2"	1	P
PRB03975	2" Silicone Port Plug	1	P
REL01408	Contactactor 2P 30A 120V Coil	1	P
REL02382	Solid State Relay 1P1T	2	P
REL03944	Blisterpak for Solid State	1	P
TRS01447	Transformer, 500VA 220V	1	P
VAL02028	Valve, CO-2 1000 PSI 120V	1	B
VAL02040	Orifice, CO-2 1000 PSI	1	P