
Job Parts List

Job # SD-608

| Part | Description | Qty Per | Pick / Buy |
|-------------|------------------------------|----------------|-------------------|
| CAS02140 | Leveling Glide Legs | 1 | P |
| CTL03235 | High Temperature Sensor | 2 | P |
| CTL03338 | PM Controller | 1 | P |
| CTL03340 | PM High/Low Limit Control | 1 | B |
| ELD03179 | Miniblock, Terminal Board | 6 | P |
| ELD03182 | Miniblock End Cap Section | 4 | P |
| ELD03286 | Quencharc Resistor, Paktron | 1 | P |
| ELE01673 | Female Connector, 9 Pin | 1 | P |
| ELF02958 | Strain Relief for 3/4" Cord | 1 | P |
| ELS01209 | Switch 2P1T | 3 | P |
| ELS01228 | Clear Cap | 1 | P |
| ELS01230 | Amber Cap | 2 | P |
| ELW02009 | Line Cord, 30A 220V | 1 | P |
| FBW01020 | Axial Fan | 1 | B |
| FBW01022 | Fan Blade, 6-1/2" SS | 1 | P |
| FBW03203 | Axial Fan Cordset | 1 | B |
| FIL02545 | CO-2 Scrubber Housing | 1 | P |
| FIL02546-1 | CO-2 Scrubber Cartridge-Each | 1 | P |
| FSE01268 | 2P 30A Fuse Block | 4 | P |
| FSE01279 | Fuse, 20A OT 250V | 2 | P |
| FSE01284 | Fuse 1A Form 101 | 1 | P |

| Part | Description | Qty Per | Pick / Buy |
|-------------|-----------------------------|----------------|-------------------|
| FSE01312 | Fuse 6-1/4A TD 250V | 1 | P |
| FSE01316 | Fuse, 25A 250V TD | 2 | P |
| FSE01351 | Fuse, 3.5A TD 250V | 2 | P |
| FSE02604 | Fuse 1A Glass | 2 | P |
| FSE02605 | Fuse Holder - Inline | 2 | P |
| GSK01162 | High Temperature Gasket | 9 | P |
| HDW01110 | Handle, Door | 1 | P |
| HDW01118 | 4 X 4 Riveted Pin Hinge | 2 | P |
| HDW01129 | Door Latch - Large | 1 | P |
| HDW01141 | Door Strike - Large | 1 | P |
| HTR01245 | Ceramic Heater - 1000W 240V | 3 | P |
| HTR03976 | Ceramic Insulators | 4 | P |
| MTR01001 | Motor 6-5/8 Shaft | 1 | P |
| PRB01582 | Port Sleeve, 1/2" | 2 | P |
| PRB01583 | Port Sleeve, 7/8" | 2 | P |
| PRB01598 | Port Sleeve 2" OD | 1 | P |
| PRB01621 | 7/8" Port Collar | 3 | P |
| PRB03975 | 2" Silicone Port Plug | 1 | P |
| REL01408 | Contactora 2P 30A 120V Coil | 1 | P |
| REL02382 | Solid State Relay 1P1T | 2 | P |
| REL03944 | Blisterpak for Solid State | 1 | P |
| TRS01447 | Transformer, 500VA 220V | 1 | P |
| VAL02028 | Valve, CO-2 1000 PSI 120V | 1 | B |
| VAL02040 | Orifice, CO-2 1000 PSI | 1 | P |

| Part | Description | Qty Per | Pick / Buy |
|-------------|-----------------------------|----------------|-------------------|
| VAL02996 | Pipe Nipple, SS 1/8" X 2" L | 1 | P |